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Beta Bond 7th generation

In the self-etch system, primer, conditioner and adhesion are done simultaneously in one solution. Due to the high amount of minerals and its unique structure, the enamel surface creates a suitable mechanical bond with this type of bonding and covers the dental cavity well.

The unique advantage of self-etching is to increase the long-term durability of the connection to dentin. This generation of bonding prevents the creation of micro-leakage due to optimal coverage and subsequently reduces sensitivity after restoration and secondary decay.

Indications:

- Suitable for all restoration classes
- To bond ceramic, veneer, crown and inlay and other indirect bonding
- Suitable for long-term restorations

Properties:

- Excellent performance in the moisture surface of dentin and suitable to use inside the canal
- Complete polymerization even in the apical area
- Contains solvents and homogeneous monomers to increase penetrability and create strong bonding
- Optimal viscosity to create minimal thickness

How to use:

The initial assessment is done to check the restoration area. The area around the damaged tooth is numbed and the tooth is ground and ready to fill.

The corners of the line and the points of ground tooth cavity should be rounded. No amalgam or other base material should remain in the prepared internal form that would interfere with light transmission and thus hardening of the restorative material.

Clean the cavity sufficiently to prevent poor bonding. If the inner surface of the cavity is contaminated with saliva or blood, wash it thoroughly and dry it completely before use of bonding.

To protect the pulp in very deep decay, you should cover around the pulp with a thin layer of the underlying filling material, calcium hydroxide, and with a suitable material (such as glass ionomer cement).

Pour a sufficient amount of bonding on the pad and place it on the inner wall of the cavity using a disposable micro brush. Remove the excess material with dry, oil-free air so that the bonding is spread as a thin layer on the surface, make sure that the entire surface of the tooth and enamel is moistened with the bonding, do not mix the bonding with other materials in any way.

| Light intensity of LED halogen lamp | Polymerization time |
|---|---------------------|
| 500mW/cm ² 1200mW/cm ² | 10 seconds |

The tip of the light cure device should be kept as close as possible and perpendicular to the surface of the resin and light for 10 seconds.

The next step is to place the composite filling material to complete the restoration process according to the manufacturer's instructions.

Maintenance instructions:

- ▶ Close the cap of the bottle immediately after each use of bonding.
- ▶ Store at 2 to 25 degrees Centigrade.
- ▶ Keep away from direct sunlight and direct light sources.
- ▶ Do not use the product after the expiration date.
- ▶ If you keep the bonding in the refrigerator, you should put it at normal room temperature for a while before using it.
- ▶ The lifespan of the product is 3 years.

Precautions:

This product is made for use in dental restorations and should be used according to the manufacturer's instructions. In case of contact with oral mucosa, rinse with water and do not use for people who are allergic to the components of this substance in any way, and if any allergic and skin sensitivities occur after using this substance, consult a specialist doctor.

Phenolic materials such as eugenol can cause a lack of proper cure, so avoid using eugenol-containing cements as liners and in the vicinity of the tooth being restored.