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
Includes: Putty - Light Body– Catalyst

Produced in accordance with the standard ISO 4823: 2015

Usage: Impression of all-ceramic bridges and crowns, orthodontic models, dental protectors during flask prosthesis

Terradent putty: With high strength and durability after hardening with very low volume changes.

Light Body: With low viscosity for accurate recording and transferring details of soft tissues and cut surfaces of teeth.

How to use: First, dry the gypsum model and the surface on which the impression will be done and prevent the flow of liquids on it.

Putty: Measure the dough precisely using the measuring cup, shape it in your palm, and press the cup into the center of the dough to create a circular groove.

Catalyst: For each cup, use an amount of catalyst equal to the diameter of the cup in length. Using too much or too little catalyst will alter the setting time of the dough.

Note: Knead the dough for 30 seconds until it becomes completely uniform in color and texture.

Impression Taking (Putty): Place the tray containing the putty into the patient's mouth and apply appropriate pressure. Wait approximately 3 minutes for the putty to set completely. (make sure the putty has fully hardened inside the patient's mouth.).

Note: High temperature and excessive amount of catalyst will shorten the setting time. Conversely, low temperature and insufficient catalyst than recommended will prolong the setting time.

Light Body: Pour an equal strip of light body and catalyst onto a glass or mixing pad, then mix thoroughly with a spatula until a homogeneous blend is achieved. Ensure that no bubbles remain during mixing. Close the catalyst lid immediately.

Impression Taking (Light Body): Place the mixture of light body and catalyst on a tray with a spatula then place it inside the patient's mouth to sets after 3 minutes.

Disinfection of the impression: Rinse the obtained impression under water, then immerse in 2% glutaraldehyde solution or sodium hypochlorite (0.5-1%) for 15 to 30 minutes.

Making Dental Cast: Terradent can be used for clinical procedures and cast preparation 30 minutes to 72 hours after impression taking. This product is compatible with dental plasters for casting.

Electroplate: Terradent is compatible with various copper and silver electroplating baths. Use only according to the electroplating system manufacturer's instructions.

Storage Conditions and Shelf Life: Store Terradent at room temperature (15 to 30 degrees Celsius) and away from direct sunlight. Under the above conditions, Terradent remains usable for up to 3 years. Always close the tube and container lid tightly after each use.


Safety Recommendations: Wash hands with soap and water after using the silicone impression material. In case of eye contact with the catalyst, rinse thoroughly with water and seek medical attention. Avoid inhaling the catalyst. Keep out of reach of children. Dispose of waste materials in designated areas.

Note: Terradent is a type of condensation silicone material. For optimal quality and performance, avoid mixing with other brands of silicone impression materials.


Technical Information	Terradent	Terradent
ISO 4823:2015	Putty consistency Type0	Light- Bodied consistency Type3
Mixing time	30 "	30 "
Working time	1' 20 "	2 '
Oral setting time	3 '	3' 30 "
Strain in compression	2-5%	2-8%
Recovery from deformation	>99%	>99%
Linear dimensional changes after 24 hours	<-0.2%	<-0.7%
Reproduction of details	20 μm	20 μm

Measurement at (23 ± 2) °c and (50 ± 10)% relative humidity

Duration of work in the initial impression

Duration of work in impression with light body

Terradent - Compression silicone for dental impression