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# Acrosun

## Self-cure acrylic to repair prosthetic restoration with Acrosun brand

### Self-cure acrylic

Acrylic powder is made of poly methyl methacrylate compound with particle size of 10 to 100  $\mu$  with other additives.

Self-cure acrylic monomer (Instant monomer) contains methyl methacrylate liquid with addition of curing agent and stabilizer.

### Product application:

- Repair of prosthesis cracks and fractures
- Reline and creating a lining on the surface of the prosthesis
- Making of partial prostheses based on metal frame

### Product features:

- Suitable working time for the usage (application) of polymer paste
- No porosity
- Lack of shrinkage
- Strong chemical bonding with the surface of the prosthesis
- Excellent polish-ability
- Proper strength

### Important notes:


- Prepare a suitable weighing tool for mixing and proceed according to the company's instructions.
- The suitable temperature for working with self-cure acrylic is 25°C in the clinic and laboratory environment.
- If the monomer is placed in the refrigerator, the setting time and working time will be longer.
- Make sure that the monomer is first poured into the mixing container and then the powder is slowly added.
- For better results, it is recommended to avoid mixing powder and monomer with other brands.

### How to prepare:


1- First, in order to create a strong bond between the self-cure acrylic and the prosthesis, we make a mechanical retainer at the repair (crack, fracture) or reline place.

2- For every 1 ml of liquid, add 2 grams of powder and pour it into the mixing bowl.

3- Mix the powder and liquid well in the mixing container for about 30 seconds and close the lid of the container.

4- After 2 minutes, the acrylic paste is ready to use and the working time for using the paste is between 3 and 4 minutes.

5- After 10 to 14 minutes, depending on the environment temperature, the polymerization is completed and the polymer paste hardens and is ready for polish.

### Self-cure monomer storage conditions:

- The liquid is flammable. Do not use near fire flames or flammable sources.
- Avoid contact with skin and eyes, in case of contact with eyes, wash immediately with water and visit a doctor.
- Do not mix with other monomers (different brands).
- Store at a temperature of  $24 \pm 4$  °C.
- Close the lid completely after each use.
- Avoid placing the liquid in direct sunlight.
- For more information on product health and safety, refer to the MSDS, which is available through the company's website.

### Powder storage conditions:

- Store in a dry and dark place.
- Keep away from heat and direct sunlight.
- It should be kept at a temperature of 15-27°C and a standard humidity of  $32 \pm 2\%$ .
- Do not mix with other powders (different brands).

Setting time	2 minutes
Working time	4 minutes
Colors available	A27- A45- A63- A81
Powder weight	500 gr
Liquid Volume	250 ml