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Acrosun

Acryl for special tray with Acrosun brand

The main base of the tray powder is poly methyl methacrylate with a particle size of 10 to 150 micrometer, which is combined with additives Silica, baking agent and pigment.

Tray monomer contains methyl methacrylate with additives such as initiator, cross-linking agent, stabilizer and anti-bubble.

Product application:

- Making special tray
- Making baseplate

Product features:


- Suitable working time for forming several bases at the same time
- Absence of porosity
- No shrinkage cage after forming on the base
- Has 2 colors, green and yellow
- Contains oil for non-adhesion to the plaster cast
- Appropriate strength

Important notes:

- Prepare accurate weighing tools and suitable containers for mixing and act according to the company's instructions. The quality of the tray is directly related to compliance with this.
- The suitable temperature to make special tray is 25 °C in a laboratory environment.
- If you put the monomer in the refrigerator, the setting time and working time will be longer.
- Be careful to first pour the liquid into the mixing container and then slowly add the powder.
- The use of powder and liquid with different brands can reduce the quality of the manufactured product.

How to prepare:


1- Apply some of the separating solution (Betafilm) carefully with a brush on the surface of the plaster cast.

2- Weigh 10 ml of tray monomer and 25 gr of powder, stir well for 30 s and cover the lid of the container.

Calculate and use the amount of materials according to the number of casts


10ml+25gr


3- After about 5 minutes, the polymer paste is ready to be removed.

Knead the paste and shape it into a 4 mm thick sheet and form it on the plaster cast.

Setting time and working time have a direct relationship with the ambient temperature, as well as the amount of powder and liquid used, if the company's instructions are followed, the application of the product will be standard.

4- After 10-15 minutes at a temperature of 25°C in the laboratory environment, the polymer paste has hardened and you can separate the special tray from the plaster cast.

Setting time	5 minutes
Working time	8 minutes
Colors available	yellow-green
Powder weight	1 kg
Liquid volume	500 ml

Tray monomer storage conditions:

- The liquid is flammable. Do not use near fire flames or flammable sources.
- Avoid contact with skin and eyes, in case of contact with eyes, wash immediately with water and see a doctor.
- Do not mix with other monomers (different brands).
- To be kept at a temperature of $24 \pm 4^{\circ}\text{C}$.
- Close the lid completely after each use.
- Avoid placing the liquid in direct sunlight.
- For more information about the health and safety of the product, refer to the MSDS, which is available through the company's website.

Powder storage conditions:

- Store in a dry and dark place.
- Keep away from heat and direct sunlight.
- Keep at a temperature of 15-27°C and a standard humidity of $32 \pm 2\%$.
- Do not mix with other powders (different brands).